


www.inventbiomed.com



Precision to Enhance Life...

INVENTO[®]-NC

RAPID EXCHANGE PTCA POST-DILATATION CATHETER



Technically Optimized for Post-dilatation

The Invento[®]-NC Rx is a High Pressure Capability and controlled expansion of non-compliant Balloon. It's a true Noncompliant Balloon with an expansion of 5%.

TECHNICAL PARAMETERS

Non Compliant polyamide Balloon	: Precise & Predictable balloon growth
Hydrophilic coating	: On the distal Shaft
Radio Opaque markers	: Zero Profile
Maximum recommended Guidewire	: 0.014" (0.36 mm)
Tip Entry Profile	: 0.016" (0.40 mm)

Excellent deflation time	: Average less than ≤ 15 sec.
Proximal Shaft diameter	: Proximal 2.19 F (0.73mm)
Distal Shaft diameter	: Distal 2.64 F (0.89 mm)
Usable Catheter length	: 1414 mm
Rapid Exchange Length	: 250 mm

ORDERING INFORMATION

Length (mm)	2.00	2.50	2.75	3.00	3.50	4.00
10 mm	INV NC 200/10	INV NC 250/10	INV NC 275/10	INV NC 300/10	INV NC 350/10	INV NC 400/10
12 mm	INV NC 200/12	INV NC 250/12	INV NC 275/12	INV NC 300/12	INV NC 350/12	INV NC 400/12
15 mm	INV NC 200/15	INV NC 250/15	INV NC 275/15	INV NC 300/15	INV NC 350/15	INV NC 400/15

Other sizes can be made available on demand.

COMPLIANCE CHART

Pressure Atm	Balloon diameter (mm)					
	2.00	2.50	2.75	3.00	3.50	4.00
6	1.85	2.32	2.56	2.77	3.24	3.67
8	1.91	2.39	2.59	2.85	3.33	3.78
10	1.96	2.44	2.69	2.91	3.40	3.87
* 12	2.01	2.50	2.74	2.98	3.48	3.98
14	2.05	2.54	2.76	3.03	3.54	4.04
16	2.08	2.57	2.79	3.08	3.58	4.11
18	2.11	2.60	2.81	3.11	3.62	4.15
** 20	2.14	2.62	2.83	3.13	3.66	4.19
22	2.16	2.65	2.85	3.16	3.68	4.23

* NP	Nominal Pressure
** RBP	Rated Burst pressure - Do not exceed